

## Endoplasmic reticulum aminopeptidase 1

Catalog No: #AP82612



Package Size: #AP82612-1 50ug #AP82612-2 100ug #AP82612-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Endoplasmic reticulum aminopeptidase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ARTS-1,Adipocyte-derived leucine aminopeptidase,A-LAP,Aminopeptidase PILS,Puromycin-insensitive leucyl-specific aminopeptidase,PILS-AP,VEGF-induced aminopeptidase,Appils,Arts1
Accession No.	Uniprot ID:Q9EQH2
Uniprot	Q9EQH2
GeneID	80898;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Erap1

Note: This product is for in vitro research use only