

Receptor tyrosine-protein kinase erbB-4

Catalog No: #AP82613



Package Size: #AP82613-1 50ug #AP82613-2 100ug #AP82613-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Receptor tyrosine-protein kinase erbB-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Proto-oncogene-like protein c-ErbB-4, Tyro-2
Accession No.	Uniprot ID:Q62956
Uniprot	Q62956
GeneID	59323;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ErbB4

Note: This product is for in vitro research use only