

Elongin-A

Catalog No: #AP82615



Package Size: #AP82615-1 50ug #AP82615-2 100ug #AP82615-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Elongin-A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	EloA,Elongin 110 kDa subunit,RNA polymerase II transcription factor SIII subunit A1,SIII p110,Transcription elongation factor B polypeptide 3,Tceb3
Accession No.	Uniprot ID:Q8CB77
Uniprot	Q8CB77
GeneID	27224;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Eloa

Note: This product is for in vitro research use only