

## Elongation of very long chain fatty acids protein 1

Catalog No: #AP82616



Package Size: #AP82616-1 50ug #AP82616-2 100ug #AP82616-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Elongation of very long chain fatty acids protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	3-keto acyl-CoA synthase ELOVL1,ELOVL fatty acid elongase 1,ELOVL FA elongase 1,Very long chain 3-ketoacyl-CoA synthase 1,Very long chain 3-oxoacyl-CoA synthase 1,SSC1,CGI-88
Accession No.	Uniprot ID:Q9BW60
Uniprot	Q9BW60
GeneID	64834;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ELOVL1

Note: This product is for in vitro research use only