

Dynamin-1-like protein

Catalog No: #AP82619

Package Size: #AP82619-1 50ug #AP82619-2 100ug #AP82619-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dynamin-1-like protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Dnm1p/Vps1p-like protein,DVLP,Dynamin family member proline-rich carboxyl-terminal domain less,Dymple,Dynamin-like protein,Dynamin-like protein 4,Dynamin-like protein IV,HdynIV,Dynamin-related protein 1,DLP1,DRP1
Accession No.	Uniprot ID:O00429
Uniprot	O00429
GeneID	10059;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DNM1L

Note: This product is for in vitro research use only