

Desert hedgehog protein

Catalog No: #AP82625



Package Size: #AP82625-1 50ug #AP82625-2 100ug #AP82625-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Desert hedgehog protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	DHH,HHG-3
Accession No.	Uniprot ID:O43323
Uniprot	O43323
GeneID	50846;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DHH

Note: This product is for in vitro research use only