

Dihydropyrimidinase-related protein 2

Catalog No: #AP82626



Package Size: #AP82626-1 50ug #AP82626-2 100ug #AP82626-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dihydropyrimidinase-related protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	DRP-2, Collapsin response mediator protein CRMP-62
Accession No.	Uniprot ID:Q90635
Uniprot	Q90635
GeneID	395155;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: DPYSL2

Note: This product is for in vitro research use only