

Dual specificity protein phosphatase 3

Catalog No: #AP82628



Package Size: #AP82628-1 50ug #AP82628-2 100ug #AP82628-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dual specificity protein phosphatase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	T-DSP11,Vaccinia H1-related phosphatase,VHR
Accession No.	Uniprot ID:Q9D7X3
Uniprot	Q9D7X3
GeneID	72349;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dusp3

Note: This product is for in vitro research use only