

DnaJ homolog subfamily C member 27

Catalog No: #AP82633



Package Size: #AP82633-1 50ug #AP82633-2 100ug #AP82633-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DnaJ homolog subfamily C member 27
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Rab and DnaJ domain-containing protein,RBJ
Accession No.	Uniprot ID:A5D7F5
Uniprot	A5D7F5
GeneID	540033;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DNAJC27

Note: This product is for in vitro research use only