

Ectonucleoside triphosphate diphosphohydrolase 5

Catalog No: #AP82637



Package Size: #AP82637-1 50ug #AP82637-2 100ug #AP82637-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ectonucleoside triphosphate diphosphohydrolase 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	NTPDase 5,CD39 antigen-like 4,ER-UDPase,Guanosine-diphosphatase ENTPD5,GDPase ENTPD5,Nucleoside diphosphatase,Uridine-diphosphatase ENTPD5,UDPase ENTPD5,Cd3914
Accession No.	Uniprot ID:Q9WUZ9
Uniprot	Q9WUZ9
GeneID	12499;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Entpd5

Note: This product is for in vitro research use only