

Histone acetyltransferase p300

Catalog No: #AP82638



Package Size: #AP82638-1 50ug #AP82638-2 100ug #AP82638-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Histone acetyltransferase p300
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	p300 HAT,E1A-associated protein p300,P300
Accession No.	Uniprot ID:Q09472
Uniprot	Q09472
GeneID	2033;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EP300

Note: This product is for in vitro research use only