

## Histone acetyltransferase p300

Catalog No: #AP82639



Package Size: #AP82639-1 50ug #AP82639-2 100ug #AP82639-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Histone acetyltransferase p300
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	p300 HAT,E1A-associated protein p300,P300
Accession No.	Uniprot ID:B2RWS6
Uniprot	B2RWS6
GeneID	328572;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ep300

Note: This product is for in vitro research use only