

## Endothelial PAS domain-containing protein 1

Catalog No: #AP82643



Package Size: #AP82643-1 50ug #AP82643-2 100ug #AP82643-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Endothelial PAS domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	EPAS-1,HIF-1-alpha-like factor,HLF,HIF-related factor,HRF,Hypoxia-inducible factor 2-alpha,HIF-2-alpha,Hif2a
Accession No.	Uniprot ID:P97481
Uniprot	P97481
GeneID	13819;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Epas1

Note: This product is for in vitro research use only