

Deoxyribonuclease gamma

Catalog No: #AP82651



Package Size: #AP82651-1 50ug #AP82651-2 100ug #AP82651-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Deoxyribonuclease gamma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	DNase gamma, DNase I homolog protein DHP2, Deoxyribonuclease I-like 3, DNase I-like 3, Liver and spleen DNase, LS-DNase, DHP2, DNAS1L3
Accession No.	Uniprot ID:Q13609
Uniprot	Q13609
GeneID	1776;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0) gene_full_name: DNASE1L3

Note: This product is for in vitro research use only