

Endothelial protein C receptor

Catalog No: #AP82654



Package Size: #AP82654-1 50ug #AP82654-2 100ug #AP82654-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Endothelial protein C receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Activated protein C receptor,APC receptor,Centrocyclin,Centrosomal protein CCD41,Endothelial cell protein C receptor,CD201,Epcr
Accession No.	Uniprot ID:Q64695
Uniprot	Q64695
GeneID	19124;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Procr

Note: This product is for in vitro research use only