

## Focal adhesion kinase 1

Catalog No: #AP82658



Package Size: #AP82658-1 50ug #AP82658-2 100ug #AP82658-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Focal adhesion kinase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	FADK 1,Focal adhesion kinase-related nonkinase,FRNK,Protein-tyrosine kinase 2,p125FAK,pp125FAK,Fak,Fak1
Accession No.	Uniprot ID:O35346
Uniprot	O35346
GeneID	25614;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PtK2

Note: This product is for in vitro research use only