

Dysbindin

Catalog No: #AP82664



Package Size: #AP82664-1 50ug #AP82664-2 100ug #AP82664-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dysbindin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Biogenesis of lysosome-related organelles complex 1 subunit 8,BLOC-1 subunit 8,Dysbindin-1,Dystrobrevin-binding protein 1,Hermansky-Pudlak syndrome 7 protein,HPS7 protein,BLOC1S8,My031
Accession No.	Uniprot ID:Q96EV8
Uniprot	Q96EV8
GeneID	84062;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DTNBP1

Note: This product is for in vitro research use only