

Developmental pluripotency-associated protein 2

Catalog No: #AP82666



Package Size: #AP82666-1 50ug #AP82666-2 100ug #AP82666-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Developmental pluripotency-associated protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Embryonic stem cell-associated transcript 15-2 protein,ECAT15-2,Pluripotent embryonic stem cell-related gene 1 protein,Phsecrg1
Accession No.	Uniprot ID:Q9CWH0
Uniprot	Q9CWH0
GeneID	73703;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dppa2

Note: This product is for in vitro research use only