

F-box only protein 38

Catalog No: #AP82668



Package Size: #AP82668-1 50ug #AP82668-2 100ug #AP82668-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	F-box only protein 38
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Modulator of KLF7 activity homolog,MoKA,SP329
Accession No.	Uniprot ID:Q6PIJ6
Uniprot	Q6PIJ6
GeneID	81545;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FBXO38

Note: This product is for in vitro research use only