

Fibroblast growth factor 5

Catalog No: #AP82670



Package Size: #AP82670-1 50ug #AP82670-2 100ug #AP82670-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fibroblast growth factor 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	FGF-5
Accession No.	Uniprot ID:Q20FD0
Uniprot	Q20FD0
GeneID	608459;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: FGF5

Note: This product is for in vitro research use only