

## Four and a half LIM domains protein 2

Catalog No: #AP82671



Package Size: #AP82671-1 50ug #AP82671-2 100ug #AP82671-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Four and a half LIM domains protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FHL-2,LIM domain protein DRAL,Skeletal muscle LIM-protein 3,SLIM-3,DRAL,SLIM3
Accession No.	Uniprot ID:Q14192
Uniprot	Q14192
GeneID	2274;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:FHL2

Note: This product is for in vitro research use only