

C-X-C motif chemokine 10

Catalog No: #AP82673



Package Size: #AP82673-1 50ug #AP82673-2 100ug #AP82673-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C motif chemokine 10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	10 kDa interferon gamma-induced protein, Gamma-IP10, Interferon-inducible protein 10, Protein Mob-1, Small-inducible cytokine B10, Inp10, Mob1, Scyb10
Accession No.	Uniprot ID:P48973
Uniprot	P48973
GeneID	245920;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cxcl10

Note: This product is for in vitro research use only