

DnaJ homolog subfamily C member 27

Catalog No: #AP82674



Package Size: #AP82674-1 50ug #AP82674-2 100ug #AP82674-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DnaJ homolog subfamily C member 27
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Rab and DnaJ domain-containing protein,Rabj,Rbj
Accession No.	Uniprot ID:Q8CFP6
Uniprot	Q8CFP6
GeneID	217378;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dnajc27

Note: This product is for in vitro research use only