

DNA polymerase alpha subunit B

Catalog No: #AP82679



Package Size: #AP82679-1 50ug #AP82679-2 100ug #AP82679-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DNA polymerase alpha subunit B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	DNA polymerase alpha 70 kDa subunit
Accession No.	Uniprot ID:Q58D13
Uniprot	Q58D13
GeneID	514793;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: POLA2

Note: This product is for in vitro research use only