

Endothelin B receptor

Catalog No: #AP82682



Package Size: #AP82682-1 50ug #AP82682-2 100ug #AP82682-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Endothelin B receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	ET-B,ET-BR,Endothelin receptor non-selective type
Accession No.	Uniprot ID:P56497
Uniprot	P56497
GeneID	403862;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EDNRB

Note: This product is for in vitro research use only