

Early growth response protein 1

Catalog No: #AP82683



Package Size: #AP82683-1 50ug #AP82683-2 100ug #AP82683-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Early growth response protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	EGR-1,Nerve growth factor-induced protein A,NGFI-A,Transcription factor Zif268,Zinc finger protein Krox-24,Egr-1,Krox-24
Accession No.	Uniprot ID:P08046
Uniprot	P08046
GeneID	13653;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Egr1

Note: This product is for in vitro research use only