

DNA polymerase alpha subunit B

Catalog No: #AP82685



Package Size: #AP82685-1 50ug #AP82685-2 100ug #AP82685-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	DNA polymerase alpha subunit B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	DNA polymerase alpha 70 kDa subunit,DNA polymerase subunit II
Accession No.	Uniprot ID:O89043
Uniprot	O89043
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pola2

Note: This product is for in vitro research use only