

S-adenosylmethionine decarboxylase proenzyme

Catalog No: #AP82688



Package Size: #AP82688-1 50ug #AP82688-2 100ug #AP82688-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	S-adenosylmethionine decarboxylase proenzyme
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	AdoMetDC,AMD
Accession No.	Uniprot ID:P17707
Uniprot	P17707
GeneID	262;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AMD1

Note: This product is for in vitro research use only