

Alpha-defensin 5

Catalog No: #AP82691



Package Size: #AP82691-1 50ug #AP82691-2 100ug #AP82691-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alpha-defensin 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Defensin-related cryptdin-5,Defcr5
Accession No.	Uniprot ID:P28312
Uniprot	P28312
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Defa5

Note: This product is for in vitro research use only