

Myotonin-protein kinase

Catalog No: #AP82692



Package Size: #AP82692-1 50ug #AP82692-2 100ug #AP82692-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myotonin-protein kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MT-PK,DM-kinase,DMK,DM1 protein kinase,DMPK,Myotonic dystrophy protein kinase,DM1PK,MDPK
Accession No.	Uniprot ID:Q09013
Uniprot	Q09013
GeneID	1760;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DMPK

Note: This product is for in vitro research use only