

C-X-C motif chemokine 9

Catalog No: #AP82695



Package Size: #AP82695-1 50ug #AP82695-2 100ug #AP82695-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	C-X-C motif chemokine 9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Accession No.	Uniprot ID:Q8K4B1
Uniprot	Q8K4B1
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CXCL9

Note: This product is for in vitro research use only