

C-X-C chemokine receptor type 7

Catalog No: #AP82696



Package Size: #AP82696-1 50ug #AP82696-2 100ug #AP82696-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | C-X-C chemokine receptor type 7 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | CXC-R7,CXCR-7,Chemokine orphan receptor 1,G-protein coupled receptor RDC1 homolog,RDC-1,Cmkor1,Rdc1 |
| Accession No. | Uniprot ID:P56485 |
| Uniprot | P56485 |
| GeneID | 12778; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cxcr7

Note: This product is for in vitro research use only