

Calcium release-activated calcium channel protein 1

Catalog No: #AP82699



Package Size: #AP82699-1 50ug #AP82699-2 100ug #AP82699-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calcium release-activated calcium channel protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Protein orai-1, Transmembrane protein 142A, Cracm1, Tmem142a
Accession No.	Uniprot ID:Q8BWG9
Uniprot	Q8BWG9
GeneID	109305;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Orai1

Note: This product is for in vitro research use only