

Dentin matrix acidic phosphoprotein 1

Catalog No: #AP82701



Package Size: #AP82701-1 50ug #AP82701-2 100ug #AP82701-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dentin matrix acidic phosphoprotein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	DMP-1,AG1,Dmp
Accession No.	Uniprot ID:O55188
Uniprot	O55188
GeneID	13406;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dmp1

Note: This product is for in vitro research use only