

Elongation factor G, mitochondrial

Catalog No: #AP82702



Package Size: #AP82702-1 50ug #AP82702-2 100ug #AP82702-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Elongation factor G, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	EF-Gmt, Elongation factor G 1, mitochondrial, mEF-G 1, Elongation factor G1, Efg, Efg1, Gfm
Accession No.	Uniprot ID:Q8K0D5
Uniprot	Q8K0D5
GeneID	28030;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Gfm1

Note: This product is for in vitro research use only