

Endothelial PAS domain-containing protein 1

Catalog No: #AP82703



Package Size: #AP82703-1 50ug #AP82703-2 100ug #AP82703-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Endothelial PAS domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	EPAS-1,Hypoxia-inducible factor 2-alpha,HIF-2-alpha,Hif2a
Accession No.	Uniprot ID:Q9JHS1
Uniprot	Q9JHS1
GeneID	29452;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Epas1

Note: This product is for in vitro research use only