

DNA excision repair protein ERCC-1

Catalog No: #AP82704



Package Size: #AP82704-1 50ug #AP82704-2 100ug #AP82704-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	DNA excision repair protein ERCC-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Ercc-1
Accession No.	Uniprot ID:P07903
Uniprot	P07903
GeneID	13870;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ercc1

Note: This product is for in vitro research use only