

# Endoplasmic reticulum-Golgi intermediate compartment protein 1

Catalog No: #AP82707

Package Size: #AP82707-1 50ug #AP82707-2 100ug #AP82707-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Endoplasmic reticulum-Golgi intermediate compartment protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ER-Golgi intermediate compartment 32 kDa protein,ERGIC-32,ERGIC32,KIAA1181,HT034
Accession No.	Uniprot ID:Q969X5
Uniprot	Q969X5
GeneID	57222;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ERGIC1

Note: This product is for in vitro research use only