

Dipeptidyl peptidase 3

Catalog No: #AP82708



Package Size: #AP82708-1 50ug #AP82708-2 100ug #AP82708-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dipeptidyl peptidase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Dipeptidyl aminopeptidase III,Dipeptidyl arylamidase III,Dipeptidyl peptidase III,DPP III,Enkephalinase B
Accession No.	Uniprot ID:Q99KK7
Uniprot	Q99KK7
GeneID	75221;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dpp3

Note: This product is for in vitro research use only