

Epithelial cell adhesion molecule

Catalog No: #AP82713



Package Size: #AP82713-1 50ug #AP82713-2 100ug #AP82713-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Epithelial cell adhesion molecule
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Ep-CAM, Epithelial glycoprotein 314, EGP314, Protein D5.7A, Tumor-associated calcium signal transducer 1, CD326, Tacstd1
Accession No.	Uniprot ID:O55159
Uniprot	O55159
GeneID	171577;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Epcam

Note: This product is for in vitro research use only