

Receptor tyrosine-protein kinase erbB-3

Catalog No: #AP82714



Package Size: #AP82714-1 50ug #AP82714-2 100ug #AP82714-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Receptor tyrosine-protein kinase erbB-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Proto-oncogene-like protein c-ErbB-3, Tyrosine kinase-type cell surface receptor HER3, HER3
Accession No.	Uniprot ID:P21860
Uniprot	P21860
GeneID	2065;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ERBB3

Note: This product is for in vitro research use only