

Estrogen receptor beta

Catalog No: #AP82716



Package Size: #AP82716-1 50ug #AP82716-2 100ug #AP82716-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Estrogen receptor beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	ER-beta,Nuclear receptor subfamily 3 group A member 2,NR3A2
Accession No.	Uniprot ID:Q9XSW2
Uniprot	Q9XSW2
GeneID	396697;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ESR2

Note: This product is for in vitro research use only