

Elongation factor 2

Catalog No: #AP82722



Package Size: #AP82722-1 50ug #AP82722-2 100ug #AP82722-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Elongation factor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	EF-2,EF2
Accession No.	Uniprot ID:P13639
Uniprot	P13639
GeneID	1938;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EEF2

Note: This product is for in vitro research use only