

Endoglin

Catalog No: #AP82723



Package Size: #AP82723-1 50ug #AP82723-2 100ug #AP82723-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Endoglin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cell surface MJ7/18 antigen,CD105,Edg
Accession No.	Uniprot ID:Q63961
Uniprot	Q63961
GeneID	13805;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Eng

Note: This product is for in vitro research use only