

Delta(3,5)-Delta(2,4)-dienoyl-CoA isomerase, mitochondrial

Catalog No: #AP82726

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP82726-1 50ug #AP82726-2 100ug #AP82726-3 1mg

Description

Product Name	Delta(3,5)-Delta(2,4)-dienoyl-CoA isomerase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Accession No.	Uniprot ID:Q13011
Uniprot	Q13011
GeneID	1891;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ECH1

Note: This product is for in vitro research use only