

# Endothelin-2

Catalog No: #AP82728



Package Size: #AP82728-1 50ug #AP82728-2 100ug #AP82728-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Endothelin-2   |
| Brief Description  | Recombinant Protein  |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE                                       |
| Species Reactivity | Rt   |
| Other Names        | ET-2,Preproendothelin-2,PPET2,Vasoactive intestinal contractor,VIC |
| Accession No.      | Uniprot ID:P23943  |
| Uniprot            | P23943   |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)             |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Edn2

Note: This product is for in vitro research use only