

Inactive dipeptidyl peptidase 10

Catalog No: #AP82729



Package Size: #AP82729-1 50ug #AP82729-2 100ug #AP82729-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Inactive dipeptidyl peptidase 10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Dipeptidyl peptidase IV-related protein 3,DPRP-3,Dipeptidyl peptidase X,DPP X,Dipeptidyl peptidase-like protein 2,DPL2,DPRP3,KIAA1492
Accession No.	Uniprot ID:Q8N608
Uniprot	Q8N608
GeneID	57628;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DPP10

Note: This product is for in vitro research use only