

Desmoglein-1

Catalog No: #AP82733



Package Size: #AP82733-1 50ug #AP82733-2 100ug #AP82733-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Desmoglein-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Desmosomal glycoprotein 1, DG1
Accession No.	Uniprot ID:Q3BDI7
Uniprot	Q3BDI7
GeneID	641355;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: DSG1

Note: This product is for in vitro research use only