

Dipeptidyl peptidase 2

Catalog No: #AP82734



Package Size: #AP82734-1 50ug #AP82734-2 100ug #AP82734-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dipeptidyl peptidase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Dipeptidyl aminopeptidase II, Dipeptidyl peptidase 7, Dipeptidyl peptidase II, DPP II, Quiescent cell proline dipeptidase, Dpp2, Qpp
Accession No.	Uniprot ID:Q9ET22
Uniprot	Q9ET22
GeneID	83768;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Dpp7

Note: This product is for in vitro research use only