

Fermitin family homolog 2

Catalog No: #AP82735



Package Size: #AP82735-1 50ug #AP82735-2 100ug #AP82735-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fermitin family homolog 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Kindlin-2,Mitogen-inducible gene 2 protein,MIG-2,Pleckstrin homology domain-containing family C member 1,PH domain-containing family C member 1,KIND2,MIG2,PLEKHC1
Accession No.	Uniprot ID:Q96AC1
Uniprot	Q96AC1
GeneID	10979;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FERMT2

Note: This product is for in vitro research use only